


Demandforce

Integration Instructions for Edge



ORTHO2

The Edge integration with Demandforce allows Edge to sync information with Demandforce including transactions, patient biographical information, appointment information, and appointment confirmations.

Setup

To enable the Demandforce integration, open the Editors from the Tools tab of the Home ribbon bar. Click System Integrations under the Integrations heading. Enter your Username and Password.

Under the Demandforce heading, click the **Enable** button.



Contact

Demandforce

Demandforce
22 Fourth Street, 12th Floor
San Francisco, CA 94103
www.demandforce.com
(415) 904-8123

Edge

Ortho2
1107 Buckeye Avenue
Ames, IA 50010
www.ortho2.com
(800) 346-4504 or (515) 233-1026