

FoxFales Doctor Time Scheduling

Integration Instructions for ViewPoint



This integration allows you to create scheduling templates in the Doctor Time Scheduling software from FoxFales and then use those templates in ViewPoint Grid Scheduler. No additional ViewPoint module is required.

Setup

As long as both FoxFales Doctor Time Scheduling and ViewPoint are properly installed, there is no additional setup necessary.


Using the Integration

Create a scheduling template in FoxFales Doctor Time Scheduling and use that software to export it as XML. Be sure to note the location the XML file is saved.

Then, in ViewPoint, follow these steps.

- 1 Be sure no one else is accessing the database and ViewPoint is closed on all other workstations.
- 2 From the ViewPoint Main Menu, click **Tools & Utilities**. Then click **DataTools**.
- 3 Click **OK** on the *DataTools Warning*. These specific instructions prevent the need to consult with our support staff before running *this* DataTool. If you have any questions or concerns, please contact Ortho2 Software Support.
- 4 On the **Destructive #2** tab, select **(131) Grid Template Import – FoxFales Dr. Time**. Click **Run**.
- 5 If you are certain you are the only one accessing your ViewPoint database, click **Continue** on the warning which appears.
- 6 Enter the location of the XML file in the **Grid Import Path** field. Click **OK**.

Your template is now part of the available templates in Grid Scheduler. To edit it further or to apply this template to specific days, refer to the instructions Chapter 8: Grid Scheduler of the *ViewPoint 7 Training & Users Guide*.



Contact

Doctor Time Scheduling
FoxFales, Inc.
PO Box 1749
Nokomis, FL 34274
www.foxfales.com
(941) 480-9440

ViewPoint
Ortho2
1107 Buckeye Avenue
Ames, IA 50010
www.ortho2.com
(800) 346-4504 or (515) 233-1026