


# GAC OrthoPlex

Integration Instructions for ViewPoint



## ORTHO2

The ViewPoint integration with OrthoPlex software from GAC allows you to start the OrthoPlex software from a ViewPoint patient record and view that patient's models. These features require the ViewPoint Third Party Integration.

## Setup

Before you begin, be sure both ViewPoint and the OrthoPlex software are installed properly.


## ViewPoint

These steps must be followed on each station where you would like to use the integration.

- 1 From the ViewPoint Main Menu, click **[Tools & Utilities]** and then **[Third Party Integration]**.

Note: If you receive a CD-Key mismatch message, contact your Ortho2 Systems Consultant.



- 2 Click **[New]**.
- 3 Select **Object Integration** for the integration type.
- 4 Enter "OrthoPlex" for the **Integration Name**.
- 5 Enter "OurOrthoPlex.OurOrthoPlexIntegration" for the **Object Name**.
- 6 Click **[Add]** and **[OK]**.



## OrthoPlex

Patients will need to be entered in the OrthoPlex database with the same patient ID they are assigned in ViewPoint. There is Patient ID field to fill out when creating an order on the OrthoPlex website.

## Using the Integration

To start OrthoPlex from a ViewPoint patient record, click the **Third Party Integration**  button in the ViewPoint toolbar and select **OrthoPlex**. OrthoPlex will open to the models of the current patient. (With this integration, the 3rd Party **Update/New** button  in the patient folder has no function.)

## Contact

**OrthoPlex**  
GAC TechnoCenter  
117 Gemini Circle, Suite 403  
Birmingham, AL 35209  
www.gacorthoplex.com  
(888) 422-2376

**ViewPoint**  
Ortho2  
1107 Buckeye Avenue  
Ames, IA 50010  
www.ortho2.com  
(800) 346-4504 or (515) 233-1026

April 18, 2014