

OrthoAnalyzer

Integration Instructions for ViewPoint



The ViewPoint integration with OrthoAnalyzer software from 3Shape allows you to start the OrthoAnalyzer software from a ViewPoint patient record and view that patient's models. These features require the ViewPoint Third Party Integration.

Setup

Before you begin, be sure both ViewPoint and the OrthoAnalyzer software are installed properly.

ViewPoint

These steps must be followed on each station where you would like to use the integration.

- 1 From the ViewPoint Main Menu, click **[Tools & Utilities]** and then **[Third Party Integration]**.
Note: If you receive a CD-Key mismatch message, contact your Ortho2 Systems Consultant.
- 2 Click **[New]**.
- 3 Select **Command Line** for the integration type.
- 4 Enter "OrthoAnalyzer" as the **Integration Name**.
- 5 For the **Filename**, use [...] to navigate to the **OrthoAnalyzer.exe** and select it.
- 6 Enter "-PatientInfoConversionType=Cmdline -PatientId=" in the **Command Line Parameters** field.
- 7 Enter "*.EMODEL" as the **Command Line Suffix**. This is case sensitive.
- 8 Check Support Biographical.
- 9 Click **[Add]** and **[OK]**.

OrthoAnalyzer

Patients will need to be entered in the OrthoAnalyzer database with the same patient ID they are assigned in ViewPoint.

Using the Integration

To start OrthoAnalyzer from a ViewPoint patient record, click the Third Party Integration button in the ViewPoint toolbar and select OrthoAnalyzer. OrthoAnalyzer will open to a list from which you select to view that patient. If the patient doesn't exist, OrthoAnalyzer will still open, but will not display a list of patients. You may then add this patient as a new patient within OrthoAnalyzer. **Note:** If the patient's ID number has a space in it, the integration will break. An enhancement is currently being worked on to fix this issue.

Contact

OrthoAnalyzer

CadBlu Inc.

37 West 28th Street, 9th Floor

New York, NY 10001

www.cadbludental.com

(212) 481-8700

ViewPoint

Ortho2

1107 Buckeye Avenue

Ames, IA 50010

www.ortho2.com

(800) 346-4504 or (515) 233-1026